

HSP90B1

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF IP FC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200 IP 1:20 - 1:50

FC 1:20 - 1:50

Calculated MW:92kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human HSP90B1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

b

Synonym:

CGP; GP96; GRP94; TRA1

Catalog #:A0879

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:7184

Isotype:IgG

Swiss Prot:P14625

Purity:Affinity purification

For research use only.

Background:

HSP90 proteins are highly conserved molecular chaperones, which normally associate with other cochaperones and play important roles in folding newly synthesized proteins or stabilizing and refolding denatured proteins after stress. HSP90B1 (GP96 or GRP94) is an endoplasmic reticulum paralogue of the cytosolic HSP90. As a major ER chaperone to mediate the UPR and a master chaperone for Toll-like receptors (TLRs), HSP90b1 chaperones peptides to MHC class I molecules of dendritic cells and other antigen-presenting cells, as well as facilitating the assembly of immunoglobulin. The protein is also involved in many other bio-processes. This antibody was generated against the C-terminal region of full-length HSP90b1.

To place an order, please [Click HERE](#).