

## ENO1

---

**Reactivity:**Human Mouse

**Tested applications:**WB IHC IF FC

**Calculated MW:**434(341)

**Observed MW:**Refer to figures

**Immunogen:**

A synthetic Peptide of human ENO1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

*To place an order, please [Click HERE](#).*

**Catalog #:**A0897

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**2023

**Isotype:**IgG

**Swiss Prot:**P06733

**Purity:**Affinity purification

For research use only.