

PLA1A

Reactivity:Human Mouse

Tested applications:WB

Calculated MW:456

Observed MW:Refer to figures

Immunogen:

A synthetic Peptide of human PLA1A

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #:A0939

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:51365

Isotype:IgG

Swiss Prot:Q53H76

Purity:Affinity purification

For research use only.