

## CTHRC1

**Reactivity:** Human

**Tested applications:** WB IHC FC

**Recommended Dilution:** WB 1:500 - 1:1000 IHC 1:50 - 1:100 FC 1:20 - 1:50

**Calculated MW:** 26kDa

**Observed MW:** Refer to Figures

**Immunogen:**

A synthetic peptide of human CTHRC1

**Storage Buffer:**

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

CTHRC1; Collagen triple helix repeat-containing protein 1; Protein NMTC1;

**Catalog #:** A0940

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 115908

**Isotype:** IgG

**Swiss Prot:** Q96CG8

**Purity:** Affinity purification

For research use only.

**Background:**

This locus encodes a protein that may play a role in the cellular response to arterial injury through involvement in vascular remodeling. Mutations at this locus have been associated with Barrett esophagus and esophageal adenocarcinoma. Alternatively spliced transcript variants have been described.

*To place an order, please [Click HERE](#).*