

CTHRC1

Reactivity:Human

Tested applications:WB IHC FC

Recommended Dilution:WB 1:500 - 1:1000 IHC 1:50 - 1:100 FC 1:20 - 1:50

Calculated MW:26kDa

Observed MW:Refer to Figures

Immunogen:

A synthetic peptide of human CTHRC1

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CTHRC1; Collagen triple helix repeat-containing protein 1; Protein NMTC1;

Catalog #:A0940

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:115908

Isotype:IgG

Swiss Prot:Q96CG8

Purity:Affinity purification

For research use only.

Background:

This locus encodes a protein that may play a role in the cellular response to arterial injury through involvement in vascular remodeling. Mutations at this locus have been associated with Barrett esophagus and esophageal adenocarcinoma. Alternatively spliced transcript variants have been described.

To place an order, please [Click HERE](#).