

RBFOX3

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW:34kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant Protein of human RBFOX3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

FLJ56884; FLJ58356; FOX3; HRNBP3; RBFOX3; NEUN; FOX-3;

Background:

NeuN, also named as FOX3 and RBFOX3, is neuronal-specific nuclear protein with MW 38-50kd(~48kd). It is one of a family of 3 mammalian Fox homologues. FOX3(Or NeuN) functions in RNA-binding protein that regulates alternative splicing events.

To place an order, please [Click HERE](#).

Catalog #:A0951

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:146713

Isotype:IgG

Swiss Prot:A6NFN3

Purity:Affinity purification

For research use only.