

PGP

Reactivity:Human

Tested applications:WB IHC FC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:100 FC1:20 - 1:50

Calculated MW:34kDa

Observed MW:Refer to figures

Immunogen:

A synthetic peptide of human PGP

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #:A10002

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:283871

Isotype:IgG

Swiss Prot:A6NDG6

Purity:Affinity purification

For research use only.