

MUC16

Reactivity: Human Mouse Rat

Tested applications: IHC

Recommended Dilution: IHC 1:50 - 1:200

Calculated MW: 2353kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human MUC16

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CA125

To place an order, please [Click HERE](#).

Catalog #: A10005

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 94025

Isotype: IgG

Swiss Prot: Q8WXI7

Purity: Affinity purification

For research use only.