

NEDD8

Reactivity:Human

Tested applications:WB IHC ICC IF FC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50- 1:200 ICC 1:50- 1:200 IF 1:20- 1:100
FC 1:20- 1:50

Calculated MW:9kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human NEDD8

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

NEDD-8

Catalog #:A10008

Antibody Type:

Monoclonal Antibody

Species:Mouse

Gene ID:4738

Isotype:IgG

Swiss Prot:Q15843

Purity:Affinity purification

For research use only.

To place an order, please [Click HERE](#).