

## RPRM

---

**Reactivity:**Mouse

**Tested applications:**WB

**Recommended Dilution:**WB 1:500 - 1:2000

**Calculated MW:**12kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human RPRM

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

*To place an order, please [Click HERE](#).*

**Catalog #:**A10013

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**56475

**Isotype:**IgG

**Swiss Prot:**Q9NS64

**Purity:**Affinity purification

For research use only.