

YAP1

Reactivity:Human

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:20 - 1:200

Calculated MW:54kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human YAP1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

fn

Synonym:

AP; YAP2; YAP65; YKI

Catalog #:A1002

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:10413

Isotype:IgG

Swiss Prot:P46937

Purity:Affinity purification

For research use only.

Background:

This gene encodes a downstream nuclear effector of the Hippo signaling pathway which is involved in development, growth, repair, and homeostasis. This gene is known to play a role in the development and progression of multiple cancers as a transcriptional regulator of this signaling pathway and may function as a potential target for cancer treatment. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Aug 2013]

To place an order, please [Click HERE](#).