

## TTL

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB

**Recommended Dilution:**WB 1:1000 - 1:4000

**Calculated MW:**43kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human TTL

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

*To place an order, please [Click HERE](#).*

**Catalog #:**A10027

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**150465

**Isotype:**IgG

**Swiss Prot:**Q8NG68

**Purity:**Affinity purification

For research use only.