

## IRF9

---

**Reactivity:**Mouse

**Tested applications:**WB IHC IF IP

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:20 - 1:50 IP 1:20 - 1:50

**Calculated MW:**44kDa

**Observed MW:**Refer to Figures

**Immunogen:**

Recombinant protein of human IRF9

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

k

**Synonym:**

RF-9; ISGF3; ISGF3G; p48

*To place an order, please [Click HERE](#).*

**Catalog #:**A1004

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**10379

**Isotype:**IgG

**Swiss Prot:**Q00978

**Purity:**Affinity purification

For research use only.