

## STX8

---

**Reactivity:**Human Mouse

**Tested applications:**WB

**Recommended Dilution:**WB 1:500 - 1:2000

**Calculated MW:**27kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human STX8

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

*To place an order, please [Click HERE](#).*

**Catalog #:**A10050

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**9482

**Isotype:**IgG

**Swiss Prot:**Q9UNK0

**Purity:**Affinity purification

For research use only.