

HLA-DRB4

Reactivity:Human Mouse

Tested applications:WB

Recommended Dilution:WB 1:1000 - 1:4000

Calculated MW:30kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human HLA-DRB4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #:A10078

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:3126

Isotype:IgG

Swiss Prot:P13762

Purity:Affinity purification

For research use only.