

GRIN3B

Reactivity:Rat

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:1000

Calculated MW:113kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human GRIN3B

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #:A10081

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:116444

Isotype:IgG

Swiss Prot:O60391

Purity:Affinity purification

For research use only.