

BAX

Reactivity:Human

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:21kDa

Observed MW:Refer to figures

Immunogen:

A synthetic Peptide of human BAX

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #:A10095

Antibody Type:

Monoclonal Antibody

Species:Mouse

Gene ID:581

Isotype:IgG

Swiss Prot:Q07812

Purity:Affinity purification

For research use only.