

## BAX

---

**Reactivity:**Human

**Tested applications:**WB

**Recommended Dilution:**WB 1:500 - 1:2000

**Calculated MW:**21kDa

**Observed MW:**Refer to figures

**Immunogen:**

A synthetic Peptide of human BAX

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

*To place an order, please [Click HERE](#).*

**Catalog #:**A10095

**Antibody Type:**

Monoclonal Antibody

**Species:**Mouse

**Gene ID:**581

**Isotype:**IgG

**Swiss Prot:**Q07812

**Purity:**Affinity purification

For research use only.