

CASP9

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW:46kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human CASP9

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #:A10097

Antibody Type:

Monoclonal Antibody

Species:Mouse

Gene ID:842

Isotype:IgG

Swiss Prot:P55211

Purity:Affinity purification

For research use only.