

## COX411

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC ICC IP FC

**Recommended Dilution:**WB 1:1000 - 1:4000 IHC 1:50 - 1:200 ICC 1:50 - 1:200 IP 1:20 - 1:50  
FC 1:20 - 1:50

**Calculated MW:**20kDa

**Observed MW:**Refer to figures

**Immunogen:**

A synthetic peptide of human COX411

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**To place an order, please [Click HERE](#).**

**Catalog #:**A10098

**Antibody Type:**

Monoclonal Antibody

**Species:**Mouse

**Gene ID:**1327

**Isotype:**IgG

**Swiss Prot:**P13073

**Purity:**Affinity purification

For research use only.