

COX411

Reactivity:Human Mouse Rat

Tested applications:WB IHC ICC IP FC

Recommended Dilution:WB 1:1000 - 1:4000 IHC 1:50 - 1:200 ICC 1:50 - 1:200 IP 1:20 - 1:50
FC 1:20 - 1:50

Calculated MW:20kDa

Observed MW:Refer to figures

Immunogen:

A synthetic peptide of human COX411

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #:A10098

Antibody Type:

Monoclonal Antibody

Species:Mouse

Gene ID:1327

Isotype:IgG

Swiss Prot:P13073

Purity:Affinity purification

For research use only.