

## KEL

---

**Reactivity:**Human

**Tested applications:**WB

**Recommended Dilution:**WB 1:200 - 1:1000

**Calculated MW:**83kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human KEL

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

*To place an order, please [Click HERE](#).*

**Catalog #:**A10116

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**3792

**Isotype:**IgG

**Swiss Prot:**P23276

**Purity:**Affinity purification

For research use only.