

CLUH

Reactivity:Human Mouse

Tested applications:WB

Recommended Dilution:WB 1:1000 - 1:4000

Calculated MW:147kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human CLUH

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #:A10140

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:23277

Isotype:IgG

Swiss Prot:O75153

Purity:Affinity purification

For research use only.