

TRPM4

Reactivity:Human

Tested applications:WB

Recommended Dilution:WB 1:200 - 1:1000

Calculated MW:134kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human TRPM4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #:A10146

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:54795

Isotype:IgG

Swiss Prot:Q8TD43

Purity:Affinity purification

For research use only.