

TAS1R3

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:1000 - 1:4000

Calculated MW:93kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human TAS1R3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #:A10157

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:83756

Isotype:IgG

Swiss Prot:Q7RTX0

Purity:Affinity purification

For research use only.