

HSP90AA1

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:1000 - 1:4000

Calculated MW:85kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human HSP90AA1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #:A10167

Antibody Type:

Monoclonal Antibody

Species:Mouse

Gene ID:3320

Isotype:IgG

Swiss Prot:P07900

Purity:Affinity purification

For research use only.