

## LRP4

---

**Reactivity:**Human

**Tested applications:**WB

**Recommended Dilution:**WB 1:500 - 1:2000

**Calculated MW:**212kDa

**Observed MW:**Refer to Figures

**Immunogen:**

Recombinant protein of human LRP4

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

CLSS; CMS17; LRP-4; LRP10; MEGF7; SOST2

**Catalog #:**A10172

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**4038

**Isotype:**IgG

**Swiss Prot:**O75096

**Purity:**Affinity purification

For research use only.

**Background:**

This gene encodes a member of the low-density lipoprotein receptor-related protein family. The encoded protein may be a regulator of Wnt signaling. Mutations in this gene are associated with Cenani-Lenz syndrome.

*To place an order, please [Click HERE](#).*