

LRP4

Reactivity:Human

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:212kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human LRP4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CLSS; CMS17; LRP-4; LRP10; MEGF7; SOST2

Catalog #:A10172

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:4038

Isotype:IgG

Swiss Prot:O75096

Purity:Affinity purification

For research use only.

Background:

This gene encodes a member of the low-density lipoprotein receptor-related protein family. The encoded protein may be a regulator of Wnt signaling. Mutations in this gene are associated with Cenani-Lenz syndrome.

To place an order, please [Click HERE](#).