

FZD5

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human FZD5

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #:A10174

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:7855

Isotype:IgG

Swiss Prot:Q13467

Purity:Affinity purification

For research use only.