

ADA

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW:41kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human ADA

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ADA1; ADA

Catalog #:A1019

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:100

Isotype:IgG

Swiss Prot:P00813

Purity:Affinity purification

For research use only.

Background:

This gene encodes an enzyme that catalyzes the hydrolysis of adenosine to inosine. Various mutations have been described for this gene and have been linked to human diseases. Deficiency in this enzyme causes a form of severe combined immunodeficiency disease (SCID), in which there is dysfunction of both B and T lymphocytes with impaired cellular immunity and decreased production of immunoglobulins, whereas elevated levels of this enzyme have been associated with congenital hemolytic anemia. [provided by RefSeq, Jul 2008]

To place an order, please [Click HERE](#).