

## FKRP

---

**Reactivity:** Human

**Tested applications:** WB

**Recommended Dilution:** WB 1:1000 - 1:3000

**Calculated MW:** 55kDa

**Observed MW:** Refer to figures

**Immunogen:**

Recombinant protein of human FKRP

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Catalog #:** A10194

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 79147

**Isotype:** IgG

**Swiss Prot:** Q9H9S5

**Purity:** Affinity purification

For research use only.

**Background:**

This gene encodes a protein which is targeted to the medial Golgi apparatus and is necessary for posttranslational modification of dystroglycan. Mutations in this gene have been associated with congenital muscular dystrophy, mental retardation, and cerebellar cysts. Several alternatively spliced transcript variants of this gene have been described, but the full-length nature of some of these variants has not been determined. [provided by RefSeq, Oct 2008]

*To place an order, please [Click HERE](#).*