

SOX17

Reactivity:Human Mouse

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50- 1:200 IF 1:20- 1:100

Calculated MW:44kDa

Observed MW:Refer to Figures

Immunogen:

A synthetic peptide of human SOX17

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #:A10198

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:64321

Isotype:IgG

Swiss Prot:Q9H6I2

Purity:Affinity purification

For research use only.