

## ZSCAN4C

---

**Reactivity:**Mouse

**Tested applications:**WB

**Recommended Dilution:**WB 1:500 - 1:2000

**Calculated MW:**57kDa

**Observed MW:**Refer to Figures

**Immunogen:**

A synthetic peptide of human ZSCAN4C

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

*To place an order, please [Click HERE](#).*

**Catalog #:**A10205

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**201516

**Isotype:**IgG

**Swiss Prot:**Q80VJ6

**Purity:**Affinity purification

For research use only.