

DIAPH2

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:2000

Calculated MW: 126kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human DIAPH2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

DIA; POF; DIA2; DRF2; POF2;

Catalog #: A10209

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1730

Isotype: IgG

Swiss Prot: O60879

Purity: Affinity purification

For research use only.

Background:

The product of this gene belongs to the diaphanous subfamily of the formin homology family of proteins. This gene may play a role in the development and normal function of the ovaries. Defects in this gene have been linked to premature ovarian failure 2. Alternatively spliced transcript variants encoding different isoforms have been identified.

To place an order, please [Click HERE](#).