

LY6E

Reactivity:Human Mouse

Tested applications:WB

Recommended Dilution:WB 1:1000 - 1:2000

Calculated MW:14kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human LY6E

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

RIGE; SCA2; RIG-E; SCA-2; TSA-1;

To place an order, please [Click HERE](#).

Catalog #:A10225

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:4061

Isotype:IgG

Swiss Prot:Q16553

Purity:Affinity purification

For research use only.