

MAN2A1

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:2000

Calculated MW: 131kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human MAN2A1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MANA2; MANII; GOLIM7; AMan II;

Catalog #: A10227

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 4124

Isotype: IgG

Swiss Prot: Q16706

Purity: Affinity purification

For research use only.

Background:

This gene encodes a glycosyl hydrolase that localizes to the Golgi and catalyzes the final hydrolytic step in the asparagine-linked oligosaccharide (N-glycan) maturation pathway. Mutations in the mouse homolog of this gene have been shown to cause a systemic autoimmune disease similar to human systemic lupus erythematosus.

To place an order, please [Click HERE](#).