

MFAP2

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 21kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human MFAP2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MAGP; MAGP1; MAGP-1;

Catalog #: A10230

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 4237

Isotype: IgG

Swiss Prot: P55001

Purity: Affinity purification

For research use only.

Background:

Microfibrillar-associated protein 2 is a major antigen of elastin-associated microfibrils and a candidate for involvement in the etiology of inherited connective tissue diseases. Four transcript variants encoding two different isoforms have been found for this gene.

To place an order, please [Click HERE](#).