

PCDH1

Reactivity:Human Mouse

Tested applications:WB

Recommended Dilution:WB 1:1000 - 1:2000

Calculated MW:115kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human PCDH1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PC42; PCDH42;

Catalog #:A10234

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:5097

Isotype:IgG

Swiss Prot:Q08174

Purity:Affinity purification

For research use only.

Background:

This gene belongs to the protocadherin subfamily within the cadherin superfamily. The encoded protein is a membrane protein found at cell-cell boundaries. It is involved in neural cell adhesion, suggesting a possible role in neuronal development. The protein includes an extracellular region, containing 7 cadherin-like domains, a transmembrane region and a C-terminal cytoplasmic region. Cells expressing the protein showed cell aggregation activity. Alternative splicing occurs in this gene.

To place an order, please [Click HERE](#).