

TPD52

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:1000 - 1:2000

Calculated MW:24kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human TPD52

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

D52; N8L; PC-1; PrLZ; hD52;

To place an order, please [Click HERE](#).

Catalog #:A10254

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:7163

Isotype:IgG

Swiss Prot:P55327

Purity:Affinity purification

For research use only.