

## ENC1

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB

**Recommended Dilution:**WB 1:1000 - 1:2000

**Calculated MW:**66kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human ENC1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

NRPB; CCL28; ENC-1; PIG10; KLHL35; KLHL37; TP53I10;

**Catalog #:**A10257

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**8507

**Isotype:**IgG

**Swiss Prot:**O14682

**Purity:**Affinity purification

For research use only.

**Background:**

This gene encodes a member of the kelch-related family of actin-binding proteins. The encoded protein plays a role in the oxidative stress response as a regulator of the transcription factor Nrf2, and expression of this gene may play a role in malignant transformation. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.

*To place an order, please [Click HERE](#).*