

EIF3B

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:1000 - 1:2000

Calculated MW:92kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human EIF3B

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PRT1; EIF3S9; EIF3-ETA; EIF3-P110; EIF3-P116;

To place an order, please [Click HERE](#).

Catalog #:A10259

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:8662

Isotype:IgG

Swiss Prot:P55884

Purity:Affinity purification

For research use only.