

## SNX17

---

**Reactivity:**Human

**Tested applications:**WB

**Recommended Dilution:**WB 1:1000 - 1:2000

**Calculated MW:**53kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human SNX17

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Catalog #:**A10274

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**9784

**Isotype:**IgG

**Swiss Prot:**Q15036

**Purity:**Affinity purification

For research use only.

**Background:**

This gene encodes a member of the sorting nexin family. Members of this family contain a phox (PX) domain, which is a phosphoinositide binding domain, and are involved in intracellular trafficking. This protein does not contain a coiled coil region, like some family members, but contains a B41 domain. This protein interacts with the cytoplasmic domain of P-selectin, and may function in the intracellular trafficking of P-selectin. Alternative splicing results in multiple transcript variants encoding different isoforms.

*To place an order, please [Click HERE](#).*