

BICD2

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:1000 - 1:2000

Calculated MW:94kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human BICD2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

SMALED2; bA526D8.1;

Catalog #:A10292

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:23299

Isotype:IgG

Swiss Prot:Q8TD16

Purity:Affinity purification

For research use only.

Background:

This gene is one of two human homologs of *Drosophila bicaudal-D* and a member of the Bicoid family. It has been implicated in dynein-mediated, minus end-directed motility along microtubules. It has also been reported to be a phosphorylation target of NIMA related kinase 8. Two alternative splice variants have been described.

To place an order, please [Click HERE](#).