

ICMT

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 32kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human ICMT

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PCMT; PPMT; PCCMT; HSTE14; MST098; MSTP098;

Catalog #: A10293

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 23463

Isotype: IgG

Swiss Prot: O60725

Purity: Affinity purification

For research use only.

Background:

This gene encodes the third of three enzymes that posttranslationally modify isoprenylated C-terminal cysteine residues in certain proteins and target those proteins to the cell membrane.

This enzyme localizes to the endoplasmic reticulum. Alternative splicing may result in other transcript variants, but the biological validity of those transcripts has not been determined.

To place an order, please [Click HERE](#).