

DDAH1

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:31kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human DDAH1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

DDAH; HEL-S-16

Catalog #:A10295

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:23576

Isotype:IgG

Swiss Prot:O94760

Purity:Affinity purification

For research use only.

Background:

This gene belongs to the dimethylarginine dimethylaminohydrolase (DDAH) gene family. The encoded enzyme plays a role in nitric oxide generation by regulating cellular concentrations of methylarginines, which in turn inhibit nitric oxide synthase activity.

To place an order, please [Click HERE](#).