

DCPS

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:2000

Calculated MW: 39kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human DCPS

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ARS; DCS1; HSL1; HINT5; HINT-5; HSPC015;

To place an order, please [Click HERE](#).

Catalog #: A10299

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 28960

Isotype: IgG

Swiss Prot: Q96C86

Purity: Affinity purification

For research use only.