

## GNG2

---

**Reactivity:** Human Mouse Rat

**Tested applications:** WB

**Recommended Dilution:** WB 1:1000 - 1:2000

**Calculated MW:** 8kDa

**Observed MW:** Refer to figures

**Immunogen:**

Recombinant protein of human GNG2

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Catalog #:** A10309

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 54331

**Isotype:** IgG

**Swiss Prot:** P59768

**Purity:** Affinity purification

For research use only.

**Background:**

This gene encodes one of the gamma subunits of a guanine nucleotide-binding protein. Such proteins are involved in signaling mechanisms across membranes. Various subunits forms heterodimers which then interact with the different signal molecules.

**To place an order, please [Click HERE](#).**