

## NLGN3

---

**Reactivity:**Human Mouse

**Tested applications:**WB

**Recommended Dilution:**WB 1:1000 - 1:2000

**Calculated MW:**94kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human NLGN3

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

HNL3;

**Catalog #:**A10310

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**54413

**Isotype:**IgG

**Swiss Prot:**Q9NZ94

**Purity:**Affinity purification

For research use only.

**Background:**

This gene encodes a member of a family of neuronal cell surface proteins. Members of this family may act as splice site-specific ligands for beta-neurexins and may be involved in the formation and remodeling of central nervous system synapses. Mutations in this gene may be associated with autism and Asperger syndrome. Multiple transcript variants encoding distinct isoforms have been identified for this gene.

*To place an order, please [Click HERE](#).*