

BANK1

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:2000

Calculated MW: 89kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human BANK1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

BANK;

Catalog #: A10312

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 55024

Isotype: IgG

Swiss Prot: Q8NDB2

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a B-cell-specific scaffold protein that functions in B-cell receptor-induced calcium mobilization from intracellular stores. This protein can also promote Lyn-mediated tyrosine phosphorylation of inositol 1,4,5-trisphosphate receptors. Polymorphisms in this gene are associated with susceptibility to systemic lupus erythematosus. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).