

BANK1

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:1000 - 1:2000

Calculated MW:89kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human BANK1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

BANK;

Catalog #:A10312

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:55024

Isotype:IgG

Swiss Prot:Q8NDB2

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene is a B-cell-specific scaffold protein that functions in B-cell receptor-induced calcium mobilization from intracellular stores. This protein can also promote Lyn-mediated tyrosine phosphorylation of inositol 1,4,5-trisphosphate receptors. Polymorphisms in this gene are associated with susceptibility to systemic lupus erythematosus. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).